




MANUFACTURERS OF RAPID DIAGNOSTIC TESTS & REAGENTS

**ONE STEP hCG IS JUST ONE SMALL STEP IN YOUR BUSY SCHEDULE
RESULTS IN AS LITTLE AS 30 SECONDS**

**HIGH QUALITY • FAST • RELIABLE • SIMPLE
CLEARED BY THE US FDA**



TO ORDER:

APT-121-25	DEVICE 25 TESTS
APT-121-50	DEVICE 50 TESTS

APT-101-25	CANISTER 25 TESTS
APT-101-50	CANISTER 50 TESTS
APT-111-25	DIPSTICK POUCHED 25 TESTS
APT-111-50	DIPSTICK POUCHED 50 TESTS

ONE STEP hCG RESULTS IN AS LITTLE AS 30 SECONDS

When lab time, quality and cost effectiveness are major concerns, AZOG hCG One-Step Pregnancy Test (Urine) becomes the wise alternative in pregnancy testing. And for good reasons:

- ▶ **AZOG hCG One-Step Pregnancy Test (Urine) is very simple to use.** No secondary reagent is required.
- ▶ It only takes one step! Add three drops of urine to the cassette or dip in urine for 5 or more seconds and results are ready and easy to read in as little as 30 seconds, depending on the concentration of hCG.
- ▶ Most important, **AZOG hCG One-Step Pregnancy Test (Urine) is extremely accurate.** For greater dependability and convenience, true 25 mIU/mL sensitivity and clear results are assured even on low hCG samples.
- ▶ **AZOG hCG One-Step Pregnancy Test (Urine) is cost effective.** In addition to saving staff time, AZOG hCG One-Step Pregnancy Test is priced dramatically lower than comparable urine hCG test kits. Furthermore, 18 months room temperature stability means no more wasted refrigerator space.

It is not surprising AZOG hCG One-Step Pregnancy Test (Urine) is today's wise alternative in laboratory hCG tests.